


DVM	
Measurement date	

Patient's name		Weight (kg)	
Breed of dog		Diagnosis	
Age (date of birth)		Limb	



STEP 1: Choose the appropriate type of orthosis for your patient.

Symbol	Measurement	Standard	Mini	Nano
Y	Distance from the middle of the carpal joint to just proximal of the lateral/palpable metacarpal heads	>5.35 cm	<5.35 cm >3.60 cm	<3.60 cm >2.75 cm
X	Minimum circumference	8.5 cm	6 cm	4.5 cm
Chosen size:		Standard	Mini	Nano



STEP 2: Measure according to chosen type.

Symbol	Measure	Name			Value [cm]
A	Circumference	at the level of the carpal joint center			
B	Width Straight line measurement. Calliper may be used, if available	at the level of the widest portion of carpal joint center, measured from the dorsal aspect			
C	Circumference	4 cm	3 cm	3 cm	
		proximal to carpal joint center			
D	Circumference* If this measurement is too proximal for your patient's limb size, place "1" for the measurement	8 cm	Not applicable	6 cm	
		proximal to carpal joint center			
E1	Length Straight line measurement	from the elbow- from the most proximal aspect of the olecranon, to the carpal joint center			
E2	Length Straight line measurement	from the carpal joint center to the ground			
F	Circumference This measurement location should not be directly covering or distal to the lateral/palpable metacarpal heads	4.5 cm	3 cm	2 cm	
		distal to carpal joint center			



*This is the measurement to help decide if your patient can have 2 proximal straps and corresponds to about +1 cm proximally, so you will want to check if this area is too proximal for comfortable range of motion in that region. Placing "1" in the WimbaAPP will let the WIMBA team know that you would like a brace with a single proximal strap.
Straight line measurement - a method where the measurement is taken as a direct point-to-point straight line, rather than following the contours of the limb.

Choose WIMBA Tokens

You can order up to 3 sets of WIMBA Tokens per order.

Blocker: Lock the orthosis at a specific angle, ROM is always 0°.

Flexion (FLX): the lower angle – maximum allowed bend.

Extension (EXT): the higher angle – maximum allowed stretch.

- ☐ Violet (FLX 75 - EXT 180)
- ☐ Blue (FLX 75 - EXT 185)
- ☐ Light Blue (FLX 75 - EXT 190)
- ☐ Beige (FLX 75 - EXT 195)
- ☐ ROM 0 blocker (180)
- ☐ Custom Token Range (FLX from 75 - EXT till 205)
FLX - EXT
- ☐ Custom Blocker (75 - 205)